

ICE-Guided LAA Closure with a Conformable Foam Implant

Vivek Y. Reddy, MD^{1,2} & Petr Neuzil, MD, PhD¹

¹Homolka Hospital, Prague; ²Helmsley Electrophysiology Center, The Mount Sinai Hospital, New York

vivek.reddy@mountsinai.org



- Conflicts of Interest:
 - *Grant support and/or Consultant:* Abbott, Biosense-Webster & Boston Scientific
 - Study funded by Conformal Medical Inc.
- This presentation discusses FDA-unapproved devices.



Disclosures

- Consultant and/or Grant support:

Abbott, Ablacon*, Acutus Medical*, Affera*, Apama Medical*, Aquqheart*, Atacor*, Autonomix*, Axon*, Backbeat*, Biosense-Webster, Biosig*, Biotronik, Boston Scientific, Cardiofocus, Cardionomics, CardioNXT/AFTx, Circa Scientific*, Corvia Medical*, Dinova-Hangzhou Nuomao Medtech*, East End Medical*, EBR, EPD*, EPIX*, EpiEP*, Eximo*, Farapulse*, Fire1, Impulse Dynamics, Intershunt*, Javelin*, Kardium*, Keystone Heart*, LuxMed*, Manual Surgical Sciences*, Medlumics*, Medtronic, Middlepeak*, Newpace*, Nuvera*, Phillips, Pulse Biosciences, Sirona Medical*, Surecor*, Thermedical, Valcare*, Vizarmed*

* I have an equity stake in these companies

• I will be discussing investigational devices without FDA or CE-Mark approval.



Conformal Prague Early Feasibility Study (EFS)

Conformal Medical's **CLAAS implant**

Single-arm, single-center
ICE feasibility study

First-in-human, n=15

STUDY PROCEDURE FLOW:

1. Conscious sedation
2. ICE in RA
3. Trans-septal
4. ICE in LA
5. Fluoroscopy-/ICE-guided placement of the LAAC device
6. TEE confirmation
7. LAAC device release



Conformal LAAC Device

Conformal Foam Matrix Cup



Two sizes to address all LAA anatomies




MOUNT SINAI
SCHOOL OF
MEDICINE



MOUNT SINAI
SCHOOL OF
MEDICINE

Conformal LAAC Device

Designed to Attempt to Minimize Thrombus



Flexible tether
eliminates cable attachment site

ePTFE cover
less thrombogenic

Compliant nitinol
endoskeleton





MOUNT SINAI
SCHOOL OF
MEDICINE



MOUNT SINAI
SCHOOL OF
MEDICINE

Allows Seal Despite Off-Axis Positioning



Conformal CLAAS Procedure



MOUNT SINAI
SCHOOL OF
MEDICINE




MOUNT SINAI
SCHOOL OF
MEDICINE

Conformal Prague-EFS ICE Study

Baseline Assessment

- 77 y/o F Hx of falls (CHA₂DS₂-VASc 5)
- TEE evaluation: 17 x 20 mm
- Implant size: regular

Patient 1901-02



MOUNT SINAI
SCHOOL OF
MEDICINE




MOUNT SINAI
SCHOOL OF
MEDICINE


Conformal Prague-EFS ICE Study

Placement & Recapture

Initial Placement (too deep)



Recapture




Patient 1901-02

Conformal Prague-EFS ICE Study


Pre-Release Assessment

Tug Test



Patient 1901-02

Angio: Seal Assessment



MOUNT SINAI
SCHOOL OF
MEDICINE





MOUNT SINAI
SCHOOL OF
MEDICINE

Conformal Prague-EFS ICE Study

ICE Evaluation




Patient 1901-09




- 74 y/o F
- Hx of epistaxis (CHA₂DS₂-VASc 4)
- TEE evaluation: 24 x 27 mm
- Implant size: regular

Conformal Prague-EFS ICE Study


Angio: Position & Leak Assessment



ICE Guided: Leak Assessment



ICE Guided: Position Assessment



Patient 1901-09



MOUNT SINAI
SCHOOL OF
MEDICINE



MOUNT SINAI
SCHOOL OF
MEDICINE

Conformal Prague-EFS ICE Study

Patient Demographics

Enrolled Patients	15
Successful implants	15 (100%)
Age	71.3 ± 10.8 yrs
Gender (Male/Female)	5 / 10 (67% / 33%)
CHA₂DS₂-VASc	4.1 ± 1.7
HAS-BLED	3.4 ± 1.4



Conformal Prague-EFS ICE Study

Results

Device procedural complications	0
LAA diameters (mean / range)	11 – 28 mm
Conformal sizes used (n)	12 regular, 2 large
Placement time (min)	55.1 ± 20.6 min
Contrast (cc)	41.5 ± 14.5
FU (45D, 6M, 1Y) as of 23-Sept-2020	15, 15, 0
Leaks (45D)	1 @ 1 mm, 1 @ 2 mm
Device related thrombus (45D)	1

As of 25Sept20



MOUNT SINAI
SCHOOL OF
MEDICINE



MOUNT SINAI
SCHOOL OF
MEDICINE

Conformal Prague-EFS ICE Study

Conclusions

The Conformal LAAC device using an ICE / Conscious Sedation Protocol:

- Is feasible
- 2 sizes address a large range of LAA anatomies
- Ability to assess seal without TEE



MOUNT SINAI
SCHOOL OF
MEDICINE



MOUNT SINAI
SCHOOL OF
MEDICINE

Conformal Next Steps

ICE Protocol Additional 5- 10 patients

- Evaluate Gen 2 Delivery System
- Validate pivotal device delivery protocol
- Complete in Q1 2021

Pivotal Trial (FDA IDE)

- RCT: Conformal LAAC vs Watchman-FLX
- 1:1, n ~ 1,400
- Primary endpoint: 1-year clinical events, seal
- Secondary endpoint: 18-month stroke & systemic embolism
- First patient: Q2 2021



Thank you



MOUNT SINAI
SCHOOL OF
MEDICINE



MOUNT SINAI
SCHOOL OF
MEDICINE